



## SuperSnap™ | High Sensitivity ATP Surface Test


SuperSnap™ High-Sensitivity Surface ATP test is the most sensitive ATP test in Hygiena’s line of ATP devices. Designed to work with all Hygiena luminometers, SuperSnap detects extremely low levels of ATP. By incorporating SuperSnap into a food safety program, risks and issues are recognized in real time so corrective action can be implemented effectively.

### Benefits:

- Patented device design ensures maximum sample collection and recovery
- More sensitive liquid stable reagent detects the lowest levels of ATP
- Perfect for meeting highest hygiene standards of food and beverage industry
- Robust reagent is resistant to sanitizer effects
- 100% recyclable plastic
- 15-month shelf life at refrigerated temperature (2 – 8 °C)
- 4-week shelf life at room temperature (21 – 25 °C)

### Sensitivity:

- 4 times more sensitive than UltraSnap™ Surface ATP Tests
- 100 times more sensitive than other systems on the market when used with EnSURE™ and EnSURE™ Touch monitoring systems
- 5 times more resistant to acid based sanitizers
- 10 times more resistant to alkali based sanitizers



Catalog No.

Description

Quantity

SUS3000


SuperSnap™ High Sensitivity ATP Test

100



Find support documents, instructional videos, and more at [www.hygiena.com](http://www.hygiena.com)

## SuperSnap Procedure



See an instructional demo at [HygienaTV](#) or [HygienaTV](#)

## Related Products

### EnSURE™ Touch

#### EnSURE™ Touch Monitoring System

EnSURE Touch is the next-generation quality monitoring system. It provides rapid and accurate sanitation verification to support your food safety program. EnSURE Touch features a 5" touch screen, wireless sync technology, and cloud-based software. EnSURE Touch is designed to adapt to your workplace, providing the data you need for complex multi-location supply chain monitoring, audit and risk management, and recall prevention.



#### AquaSnap™ Total and AquaSnap™ Free

AquaSnap measures ATP found in liquid samples, utilizing a honeycomb dipper to consistently collect 100 µl of liquid sample. AquaSnap is available in Total and Free formats for monitoring microbial and non-microbial activity in water samples.



Catalog No.	Description	Quantity
ETOUCH	EnSURE™ Touch Monitoring System	1
AQ-100X	AquaSnap™ Total Water ATP Test	100
AQ-100FX	AquaSnap™ Free Water ATP Test	100



Find support documents, instructional videos, and more at [www.hygiena.com](http://www.hygiena.com)